

3.5.P07**Experimental diamond growth in carbonate-KCl: An analogue for the growth of diamond coat in the mantle**E.L. TOMLINSON AND A.P. JONES

University College London, Gower Street, London, WC1E 6BT, U.K. (emma.tomlinson@ucl.ac.uk)

High-pressure, high-temperature diamond growth experiments have been conducted in the system C-K₂CO₃-KCl at 1000-1400°C, 5.3 GPa. KCl is of interest because of the strong effect that halogens have on the phase relations of carbonate-rich systems and because of the occurrence of KCl coexisting with alkali silicate-carbonate fluids in natural coated diamond. Coated diamonds comprise a core of gem quality diamond, surrounded by diamond rich in fluid inclusions. Micro-inclusions in coated diamond samples from the Democratic Republic of Congo contain fluids that vary between three end-member compositions: K-Cl rich, Si-Al-K rich and Ca-P rich. Inclusions close to the boundary between the transparent core and the fibrous coat are particularly enriched in the K-Cl end-member, suggesting a link between KCl and the onset of coat growth. We have used systems C-K₂CO₃-KCl and C-K₂CO₃-H₂O as analogues for carbonate-KCl mantle fluids in diamond growth experiments. Diamond growth is spontaneous, and the run products are analysed by X-Ray Diffraction and Scanning Electron Microscopy.

The presence of KCl reduces liquidus temperature of potassium carbonate, allowing it to act as a solvent catalyst for diamond growth at temperatures below the continental geotherm. It is suggested that this is due to disruption of the carbonate lattice by Cl. Chlorine does not effect the diffusion rate of carbon in the carbonate melt. Diamond growth using carbonate solvent catalysts is characterised by a relatively long induction period. However, the addition of KCl also reduced the period for diamond growth in carbonate to <<5 minutes; therefore such induction period appears to be necessary.

3.5.P08**Helium-Neon systematics in OIB and the nature of the source of mantle plumes**M. MOREIRA, S. DOUCET, P.M. MADUREIRA, A. LECOMTE AND C.J. ALLÈGREInstitut de physique du globe de Paris, France
(moreira@ipgp.jussieu.fr)

We have designed a new extraction and purification line, fully automated, coupled with our glass mass spectrometer ARESIBO I, equipped with a faraday cup and an ion counting system. The very low CO₂ and ⁴⁰Ar levels of the machine allow us to analyze the neon isotopic composition in very small samples. We have performed helium and neon analyses on olivine phenocrysts from different oceanic island basalts (St Helena, Gough, Tristan da Cunha, Azores, Kerguelen, Hawaii). Most of the samples present anomalies compared to air for neon (with ratio ²⁰Ne/²²Ne up to 11.67) and all the data plot between the MORB line (e.g. Sarda et al. 1988) and the Iceland line (e.g. Moreira et al., 2001). Corrected for air ²¹Ne/²²Ne ratios (assuming a solar ²⁰Ne/²²Ne) are between 0.036 and 0.081. In a diagram ²¹Ne/²²Ne_c versus 4He/3He, all the OIB plot between the mean MORB composition and a solar like composition, and can be interpreted as a binary mixing with a parameter ($= (^3\text{He}/^{22}\text{Ne})_{\text{MORB}} / (^3\text{He}/^{22}\text{Ne})_{\text{solar}}$) for the hyperbolae around 10 suggesting that mixing between the plume material and the MORB material occur between liquids. Another interpretation is that the ³He/²²Ne of the MORB source is higher than the solar ³He/²²Ne ratios which appears in contradiction with previous estimates. Our conclusion is that the neon is a more sensitive tracer of primordial material than helium.